



Portable Conductivity/Resistivity/TDS/Salinity Meter

CDH-287-KIT



- ✓ 6 Conductivity Ranges
- ✓ 2 Resistivity Ranges
- ✓ 5 TDS Ranges
- ✓ Practical Salinity in the Range of 2 to 42.0, in Accordance with UNESCO Data
- ✓ Programmable Temperature Coefficient
- ✓ Microprocessor-Based
- ✓ Easy to Operate
- ✓ Ideal for Most Water Applications

The CDH-287 is a portable, multi-ranging conductivity meter with unsurpassed accuracy and reliability in the field. It also measures resistivity, total dissolved solids and practical salinity, making the CDH-287 the most versatile meter on the market. It comes with a glass, dip-style conductivity probe with an integral temperature sensor and cell constant of 1.0. This microprocessor-based conductivity meter features auto-ranging, programmable temperature coefficient and error diagnostics. Simply pressing the keypad switches from conductivity to TDS, resistivity or salinity units. A concentration mode allows operators to choose their own concentration units and create a specific calibration curve, by measuring up to four concentration standards and inputting the values.

To increase range accuracy, cells are also available with constants of K=0.1 and K=10. Dip cells and flow cells are available in both glass and epoxy body styles.

Specifications

Conductivity

Ranges:

0.00 to 19.99 $\mu\text{S/cm}$ or 1.999 mS/cm
 0.00 to 199.9 $\mu\text{S/cm}$ or 19.99 mS/cm
 0000 to 1999 $\mu\text{S/cm}$ or 199.9 mS/cm
 0.00 to 19.99 mS/cm or 1999 mS/m
 00.0 to 199.9 mS/cm or 19.99 mS/m
 000 to 1999 mS/cm or 199.9 mS/m

Accuracy: $\pm 0.3\%$ rdg

Temperature Compensation:

Automatic, 0 to 50°C (32 to 212°F)

Total Dissolved Solids

Ranges:

0 to 19.99 mg/L , 0 to 199.9 mg/L
 0 to 1999 mg/L , 0 to 1.999 g/L
 0 to 19.99 g/L , 0 to 132.0 g/L

Accuracy: $\pm 0.3\%$ rdg

Resistivity

Ranges:

0 to 1.999 $\text{M}\Omega/\text{cm}$ or 0.019 $\text{M}\Omega/\text{m}$
 0 to 19.99 $\text{M}\Omega/\text{cm}$ or 0.199 $\text{M}\Omega/\text{m}$

Accuracy: $\pm 0.3\%$ of reading

Salinity

Range: 2.0 to 42.0% salinity; automatic conversion from conductivity, using temperature relationship for seawater, in accordance with UNESCO, IASPO data

Accuracy: $\pm 0.3\%$ rdg

Temperature Compensation: Automatic, -2.0 to 35°C (28.4 to 95°F)

Concentration


Range: 0 to 9999, automatic ranging, choice of units, background offset function

Calibration: 4 point straight line interpolation

Temperature

Range: -30.0 to 130.0°C (-22.0 to 266.0°F)

Accuracy: $\pm 0.3^\circ\text{C}$ (0.5°F)



CDH-287 shown with the CDE-5001-GDI glass dip probe K=1.0 with ATC sensor.

Shown smaller than actual size.

General Specifications

Reference Temperature: 25°C (77°F), selectable to 20°C (68°F)

Temperature Coefficient: Preset to 2%/°C; programmable from 0 to 5%/°C

Measurement Frequency: 3000 Hz

Recorder Output: ± 200.0 mV

RS232

Display: 12.7 mm (0.5") LCD

Dimensions: 100 L x 180 W x 44 mm D (3.9 x 7.1 x 1.7")

Cable: 1.2 m (3.9')

Weight: 410 g (0.9 lb)

Power: 9V battery (included)



Shown smaller
than actual size.

CDH-287 shown with the
CDE-5001-GDI glass dip
probe K = 1.0 with
ATC sensor.

CDH-287-KIT is supplied with meter, glass
dip-style conductivity probe with integral
temperature sensor and cell constant of 1.0,
sample bottle, calibration solution, rugged
carrying case and operator's manual.

To Order

Model No.	Description
CDH-287-KIT	Conductivity/resistivity/TDS/salinity meter, glass dip-style conductivity probe with integral temperature sensor, K = 1.0, carrying case, sample bottle, 9V battery, calibration solution
CDE-5001-GDI	Replacement conductivity probe, glass dip-style, K=1.0 with ATC, platinum plates, 12 x 130 mm (0.47 x 5.1")
CDE-5002-PD1	Polymer dip-style probe, K = 1.0 with ATC, platinum plates, 12 x 130 mm (0.47 x 5.1")
CDE-5004-ED10	Epoxy dip-style probe, K = 10 with ATC, carbon plates, 26 x 353 mm (1.0 x 13.9")
CDE-5005-GF1	Glass flow cell, K = 1.0 with ATC, platinum plates, 13 x 166 mm (0.5 x 6.5") overall length, 5 mm (0.2") tubing connections, 33 mm (1.3") cell head, 4 mL minimum volume
CDE-5008-EF10	Epoxy flow cell, K = 10 with ATC, carbon plates, 26 x 203 mm (1.0 x 8.0") with 10.5 mm (0.4") tubing connections, 353 mm (13.9") overall length, 13 mL volume
CDE-5010-ED1	Epoxy dip-style probe, K = 1 with ATC, carbon plates, 26 x 250 mm (1.0 x 9.8")
CDE-5011-ED01	Epoxy dip-style probe, K = 0.1 with ATC, carbon plates, 26 x 216 mm (1.0 x 8.5")
CDE-5012-EF1	Epoxy flow cell, K = 1.0 with ATC, carbon plates, 26 x 100 mm (1.0 x 3.9") with 10.5 mm (0.4") tubing connections, 250 mm (9.8") overall length, 26 mL volume
CDE-5013-EF01	Epoxy flow cell, K = 0.1 with ATC, carbon plates, 26 x 66 mm (1.0 x 2.6") with 10.5 mm (0.4") tubing connections, 216 mm (8.5") overall length, 12 mL volume
CDE-5014-GD01	Glass dip-style probe, K = 0.1 with ATC, platinum plates, 20 x 130 mm (0.8 x 5.1")
CDE-5019-ED1	Epoxy dip-style probe, K = 1.0 with ATC, carbon plates, 12 x 110 mm (0.5 x 4.3")
MN1604	3 Replacement 9V battery

CDH-287-KIT is supplied with meter, glass dip-style conductivity probe with integral temperature sensor and cell constant of 1.0, sample bottle, 9V battery, calibration solution, rugged carrying case and operator's manual.

Ordering Examples: CDH-287-KIT, meter and accessories, plus CDE-5008-EF10, epoxy flow cell.

CDH-287-KIT, meter and accessories, plus CDE-5011-ED01, 0.1 cell constant probe.

Cell Constant	Measuring Range	Typical Applications
0.1	>100 μ S (platinum) >200 μ S (carbon)	Pure demineralized, distilled or boiler-fed water
1.0	100 μ S to 100 mS	Surface or wastewater-diluted salt solutions, fertilizers, electroplating rinses
10	Over 100 mS	Concentrated salt solutions, sea water

Note: Carbon (Graphite) probes are easier to clean and are recommended when suspended solids are present.

Order Online

omega.com

Over 100,000 Products Available!

To download information and to order Conductivity
Instrument products online, visit **omega.com**

E-16